




TECHNOLOGY MATRIX HIGH-FLOW NASAL CANNULA

Company	FLEXICARE	FLEXICARE	FLEXICARE
	 <p>NioFlo™ Nasal Cannula</p>	 <p>Veoflo™ High Flow Nasal Cannula</p>	 <p>High-Flow Nasal Cannulas</p>
Website	www.flexicare.com	www.flexicare.com	www.flexicare.com
Patient population	Pediatric, neonatal	Adult	Pediatric, neonatal
Sizes available	Micro Premature, Premature, Neonate, Infant, Pediatric Small, Pediatric Medium, Pediatric Large	Small (4mm) Medium (5mm) Large (6mm)	Premature through adult
O ₂ flow rates supported	Premature, Neonate, Infant, & Intermediate Infant: 8 LPM Pediatric: 20 LPM	Up to 60 LPM	Premature, Neonate, Infant, & Intermediate Infant: 8 LPM Pediatric: 20 LPM Pediatric/Small Adult, Adult: 40 LPM
Compatibility	Included adapter provides a secure connection to both 15mm and standard O ₂ tubing	Connectors available as an individual item or with a Veoflo® to suit your set up requirements	The included adapter provides a secure connection to both 22mm and 15mm respiratory limbs
Construction materials	Satin finished star lumen tubing, adjustable slider	Secure and Adjustable Tube Holder reduces drag and works in tandem with the lanyard and crocodile clip to provide extra support. Easy to adjust wide elastic split strap provides a secure fit and increased patient comfort. Quick and Secure Setup Wings on the swivel connector allow easy connection and disconnection. Detachable smoothbore tube can be comfortably placed on either side.	Adjustable Slider, Soft Touch Star Lumen tubing
Comfort fitting solutions	Soft, premium and satin finished materials with rounded edges that are gentle on the face; minimizing friction and pressure on delicate tissues.	Soft contoured nasal prongs aid patient compliance	Included sizing guide assists in selection. Soft Touch Tubing is gentle on the ears and face, minimizing friction and pressure on delicate tissues.
Quantity sold per package	Box of 10	Box of 25	Box of 25
Additional features	NioFlo™ nasal cannulas offer many advantages when used in conjunction with humidification, including accurate delivery of high and low flow oxygen and optimized mucociliary clearance, while increasing alveoli recruitment with a premium cannula that offers improved patient comfort and compliance. Comfortable cannula for optimal care.	Veoflo™ High-Flow Nasal Cannula is designed to allow for a range of oxygen concentrations using variable flow rates at optimal conditions. The Veoflo™ range includes a Tracheostomy Interface that is ideal for delivering high oxygen concentrations to spontaneously breathing tracheostomized patients during weaning from controlled ventilation.	Flexicare's comprehensive range of High-Flow Nasal Cannulas (HFNC) offer many advantages when used in conjunction with humidification, including accurate delivery of high FiO ₂ and optimized mucociliary clearance, while increasing alveoli recruitment in a technique that offers improved patient comfort and compliance.